



## Revolver Syringe



Weiyou® Revolver syringe is made of high quality metal. Weiyou® Revolver syringe is often used to administer injections, insert intravenous drugs into the bloodstream, apply compounds such as glue or lubricant, and measure liquids. As the manufacturer in China, Weiyou® has mature production lines and supply chains, and always supply best quality product, best favourable price and best service.

## Product Description

Weiyou® Revolver syringe is made of high quality metal. Weiyou® Revolver syringe is often used to administer injections, insert intravenous drugs into the bloodstream, apply compounds such as glue or lubricant, and measure liquids. Weiyou® Revolver syringe is mainly used for Cattle, Horse Sheep, Goat and other animals. OEM customization is acceptable. Vaccination to prevent disease for all kinds of animals. Weiyou® Revolver syringe is continuous injection and accurate scale. Free samples, and small quantity can be supplied.

## Specification

Item	Revolver syringe
Material	Metal
MOQ	100pcs
Specification	5ml, 10ml, 20ml, 30ml, 50ml
Use For	Cattle, Horse Sheep, Goat, etc
Advantage	Continuous injection, accurate scale
Type	Vaccination to prevent disease
Packing	Carton

## Features

All metal parts are made of brass material, nickel or chrome plated  
Marks on glass barrel



[www.nbweiyou.com](http://www.nbweiyou.com)

Perfect dose adjustment

Suitable for the inoculation of livestock

Clean well and oil with silicone oil, regularly replace washers



































